



CHARLOTTESVILLE-ALBEMARLE
RESCUE SQUAD



Glucose Monitor Weekly Control Documentary

(3 weeks per sheet)

Bayer Ascensia Contour Blood Glucose Monitoring System

| | | | | | |
|--------------------------------------|---|----------------------|--------------------------------------|------------------------------|---|
| Date of Test | CARS Device # | BGMS Serial # | Test Strip Lot # | Strip Expiration Date | Normal Control Range |
| <i>Ex: 10/9/07</i> | <i>Ex: #1</i> | <i>Ex: 3078148</i> | <i>Ex: 7AC3C53</i> | <i>Ex: 01/09</i> | <i>Ex: 100-138</i> |
| | | | | | |
| Normal Control Solution Lot # | Control Solution Expiration Date | Test Result | Pass or Fail? If Fail Repeat. | Repeated Test Result | Fail 2nd time? YES, Take Out of Service |
| <i>Ex: 1508</i> | <i>Ex: 03/09</i> | <i>Ex: 107</i> | <i>Ex: Pass</i> | <i>Ex: N/A</i> | <i>Ex: N/A</i> |
| | | | | | |
| XXXXX | XXXXX | XXXXX | XXXXX | XXXXX | XXXXX |
| Date of Test | CARS Device # | BGMS Serial # | Test Strip Lot # | Strip Expiration Date | Normal Control Range |
| | | | | | |
| Normal Control Solution Lot # | Control Solution Expiration Date | Test Result | Pass or Fail? If Fail Repeat. | Repeated Test Result | Fail 2nd time? YES, Take Out of Service |
| | | | | | |
| XXXXX | XXXXX | XXXXX | XXXXX | XXXXX | XXXXX |
| Date of Test | CARS Device # | BGMS Serial # | Test Strip Lot # | Strip Expiration Date | Normal Control Range |
| | | | | | |
| Normal Control Solution Lot # | Control Solution Expiration Date | Test Result | Pass or Fail? If Fail Repeat. | Repeated Test Result | Fail 2nd time? YES, Take Out of Service |
| | | | | | |